



HIP

PRODUCT BROCHURE OF HIP PROSTHESIS

HIP Replacement Solutions

COMPANY INTRODUCTION

Just Medical Devices (Tianjin) Co., Ltd., established in 1958, is a leading Chinese manufacturer of hip and knee joint implants and instruments dedicated to providing high-quality solutions for joint diseases to patients worldwide. We are integrating R&D, manufacturing, sales, training, and service to provide customers with safe, effective, and exquisite medical products and services to the fullest extent possible.

Just Medical offers six product lines: Hip Reconstruction, Knee Reconstruction, Partial Reconstruction, New Interface Prosthesis, Regenerative Repair, and Al Digitalization.

Our products have a wide range of specifications, are easy to operate, and meet the needs of various clinical procedures, including minimally invasive, primary, complex primary, primary revision, and complex revision surgeries. We provide end-to-end process, digital, stepwise product solutions for treating joint diseases.

Our products are certified by GMP, ISO, CE, and NMPA. We have obtained 33 global Class III medical device registrations and applied for 245 patents, including over 100 invention patents. With a sales network spanning 49 countries and regions globally and over 550,000 implants worldwide, we have gained the trust of doctors and patients worldwide, leading to an increasing reputation in the international market.

We possess multiple advanced product core technologies and have introduced groundbreaking innovations to the market, including the world's first clinically validated 3D-printed trabecular modular femoral stem and the world's first 3D-printed zonal trabecular knee prosthesis system, addressing market gaps. In the future, Just Medical will continue to invest in innovative platforms such as AI, 3D printing, new material applications, and regenerative medicine. Adheres to the principles of "Full product line, On-time supply, Value for money," focusing on joint step therapy to elevate the quality of human life.



CONTENT

HIP PRODUCT FAMILY	1
QUALITY CONTROL	3
Germany Raw Material ·····	3
Precise Processing ·····	3
Strict Inspection	3
Wear Test in EndoLab	4
FEMORAL STEMS	5
MINI™ Minimally Invasive Stem ······	5
DELTA™ Classic Rectangular Stem ······	6
DELTA™ Rectangular Stem ·····	7
HARMONY™ Tapered Stem ····································	8
TAICH® Cemented Stem ·····	9
TAICH® Long Cemented Revision Stem ······	10
SEE [®] Trabecular Modular Stem ······	11
RSL® Cementless Revision Stem	12
ACETABULAR CUP	13
Revision Acetabular Components ·····	13
3D SEE® Trabecular Titanium Cup	14
Highly Crosslinked Polyethylene Liner	16
HARMONY™ Acetabular Cup ······	17
HARMONY™ Liner/Constrained Liner ······	18
Dual Mobility Acetabular System	19
SAFETY [™] Bipolar Head ····································	20
Cemented Acetabular Cup ······	20
FEMORAL HEAD	21
CoCrMo/SS Head ·····	21
Ceramic Head	21
Mysmond [™] Zirconium-niobium Femoral Head ·······	21
CUSTOMIZED PROCESSES	22
HIP INSTRUMENTATION	23



HIP PRODUCTS FOR ALL-ROUND SURGICAL SOLUTIONS



Dynamic fatigue tests of femoral stem's head-neck conjunction and body after 10 million cycles in the international CNAS laboratory shows excellent results and no risk of fracture. Dynamic wear tests after 5 million cycles in the international Endolab® laboratory in Germany shows excellent wear resistance.

HARMONY™ ACETABULAR CUP SYSTEM

SEE® 3D PRINTING TRABECULAR ACETABULAR CUP SYSTEM

CAGE AND MESH SYSTEM

CEMENTED ACETABULAR CUP

ACETABULAR CUP













SEE Trabecular Cup









AOS Ring (Revision)

AOS Cage (Revision)





LINER













32/10° (HPE)



Dual Mobility Metal Liner



(46-58) Dual Mobility PE Liner



36-54 Ceramic Lining





32 Cemented Liners (HPE)



FEMORAL HEAD



Φ 22 (0/+3.5)



Φ24 (+0/3.5/7)



Φ28 (0/±3.5/+7)



Φ32 (0/±3.5/+7)

PRIMARY-



Φ28 (S/M/L)





Φ36 (S/M/L/XL)





Φ28 (S/M/L) BIOLOX[®]Delta Ceramic



Φ32 (S/M/L/XL) BIOLOX [®] Delta Ceramic



Φ36 (S/M/L/XL) BIOLOX[®]Delta Ceramio



Customized Prosthesis Design

FEMORAL STEM























-customized-REVISION -

-MINIMAL INVASIVE-

--- COMPLEX PRIMARY ----

Germany Raw Material

All raw material of UHMWPE inserts were manufactured in Germany, meeting the technical requirements in ISO 5834 part 2 and ASTM F648.



Precise Processing









Strict Inspection

JUST MEDICAL Inspection Center









Wear Test in EndoLab®

EndoLab® GmbH, a spin-off from the Technical University of Munich, Germany, offers a variety of technological implant testing services to develop and certify medical products. It is an accredited (DAkkS O-PL-18838-02-00) and certified (ZLG-P-944.98.07) test laboratory according to DIN EN ISO/IEC 17025 and 93/42/ EWG.



▲ Kinematics and load profile of the ISO 14242-1 EndoLab hip simulator.



▲ Statistical data of the mean wear rate values found for all comparative THR systems tested at EndoLab® so far (n=21). The value established for the cup liner of the JUST MED THR system is indicated by a red dot.

The UHMWPE cup liners showed a mean wear rate of 12.53 mg per million cycles (SD 1.47 mg per million cycles) after 5 million cycles. The wear rate was determined between 0 and 5 million cycles.

To date, EndoLab® has tested n=21 comparative THR systems with a CoCrMo femoral head articulating against a conventional UHMWPE Iliner (not-aged). A mean wear rate of 27.49 mg per million cycles (SD 16.03 mg per million cycles) was found. The lowest wear rate measured was 2.71 mg per million cycles, the highest wear rate was 59.69 mg per million cycles and the median was 27.16 mg per million cycles.

The mean UHMWPE cup liner wear rate of the JUST MED THR system is below the mean value of all comparative THR systems tested at EndoLab®, so far.

MINI

Minimally Invasive Stem

- Minimally invasive: soft tissue preservation and bone conservation
- Stability: initial fixation and axial & rotational stability
- Specification: 6-15#, 1mm increment in M/L and 2mm in length
- Ti+HA proximal coating provides extensive bone ingrowth



CLASSIC Rectangular Stem

- A Rectangular design for excellent torsional stability
- Reliable initial stability
- ▲ Long-term biological fixation
- Simple instrumentation

The reduced lateral shoulder eases stem insertion

Rectangular cross section design resists axial/ torsional stresses for initial stability

The proximal steps designed to convert hoop stresses to compression loads to avoid intra-operative fracture, and provide consistent implant seating height and initial stability

Distal vertical/horizontal grooves provide rotational and axial stability

Optimal offset aiming to restore hip biomechanics




Ti6Al4V material and Ti+HA double layered coating (225 micron Ti coating plus 70 micron HA coating)

ensure long-term bone ingrowth



Rectangular Stem


- Rectangular Dimension conserves bone while providing excellent torsional stability
- HA+Ti Porous Coating allows for excellent biological fixation
- Simple Instrumentation





Tapered Stem

- HA+Ti Porous Coating allows for excellent biological fixation to bone
- Reduced Neck Geometry allows for greater range of motion
- Excellent biomechanics performance








TANCH Long Cemented Revision Stem

The Threaded Hole is convenient to control rotation during implantation Polished Narrowed Neck minimizes the risk of dislocation and increases the range of motion Markings on the proximal stem 132 degrees neck-shaft is to assist positioning to restore angle exact leg length The mirror-polished stem creates radial compressive loading to reduce shear force as stem migrates distally Distal centralizer prevents the tilt of the prosthesis in the medullary cavity and equalizes thickness of the cement mantle








Trabecular Modular Stem




^{**} The proliferation and growth experiment of osteoblasts shows the 3D printing trabecular structure of titanium alloy has good biocompatibility. The microporous structure is conducive to the growth of bone cell.





Product Features

- 2° taper design provides a more even load distribution on the stem, avoiding proximal stress shielding.
- 8 longitudinal ridges offer excellent rotational stability and space for bone tissue regeneration.
- The circular cross-section design allows for intraoperative adjustment of anteversion, providing greater flexibility.
- Surface roughening of bio-compatible materials ensures excellent longterm stability.
- Proximal binding hole design accommodates patients with greater trochanter fractures.





Revision Acetabular Components



— AOS Cage —

- ▲ The groove design on the acetabular rim realizes better attachment between the prostheses and acetabular rim
- ▲ The post-lateral region of the cement buttress, provides more support for cement
- ▲ Three flanges enhance the cup fixation and provide more screw insertion choices

- ▲ Various screw holes meet requirements for different clinical cases
- ▲ Perfect rough surface for cement-bone interface





— Acetabular Mesh —



Trabecular Titanium Cup

- ▲ Titanium trabecular structure for greater biocompatibility and superior stability
- -Friction coefficient of 1.08 provides excellent initial stability
- ▲ Reliable long-term bone ingrowth
- -3D through-hole, 75%-80% porosity and 600-800µm pore size facilitate bone ingrowth
- -Modulus of elasticity similar to cancellous bone avoids stress shielding and osteolysis
- -Electron beam melting technology without risk of shedding
- ▲ Excellent Prosthesis Strength
- ▲ A variety of sizes available
- -DDH ,primary and revision type



SEE (DDH)

SEE (Primary)





▲ A specialized acetabular cup with a ceramic lining, featuring a 3D printed bone trabecular structure for precise replication of physiological bone trabecular micropores, resulting in optimal initial stability and long-term performance stability.





Trabecular Titanium Cup

- The 6 screw holes at the edge and the 6 screw holes at the dome were designed with uniform distribution, The operation was flexible, It can make the doctor choose the position of the screw more reasonable during the operation
 - Titanium trabecular bone structure has good biocompatibility
 - Excellent initial stability: 1.08 Friction coefficient, Compact pressing design, The 180° hemisphere design
 - Excellent osseointegration performance: 600-800µm pore size, 75%-80% porosity, Three dimensional through porous structure
 - The cup and liner were fixed by a card lock connection, HPE-10° liner and Harmony10° liner were used to meet the needs of acetabular revision surgery
 - Occlusion at a position without acetabular screw insertion increases the acetabular contact area and increases the overall strength of the cup



3D JCT (Revision)

▲ "Starry sky" screw hole design

- -Evenly distributed 9 screw holes in the outer edge and 6 in the dome, allowing surgeons for desired screw locations during the surgery
- -The "starry sky" design is more suitable for hip revision patients with severe acetabular defect and irregular acetabulum





Highly Crosslinked Polyethylene Liner



- ▲ German highly cross-linked polyethylene material, 10 years clinical experience
- ▲ Gama irradiation with heat treatment post-process, achieve a perfect balance in wear resistance, oxidation resistance and mechanical properties
- ▲ The clinical proven design of HARMONY 10°liner
- Available in various options, meet different patient's need

HPE Liner (10°)

Cemented Liner-HPE

▲ Design of a combination of SEE cup (Revision) + Cemented liner





Acetabular Cup

- ▲ Multiple specifications meet the need of different patients
- ▲ 0.375mm coating thickness provides the optimal press-fit fixation effect, with over 15 years proven clinical experience
- ▲ Mirror polished inner surface reduces wear and debris



-Primary-







- ▲ 180° femoral head coverage for greater range of motion
- ▲ Taper locking between the liner and cup for better transmission of force while reducing wear
- ▲ German polyethylene material reduces postoperative wear and debris

— HARMONY Liner 10° —





Dual Mobility Acetabular System



▲ Metal Liner

- CoCrMo alloy
- Mirror polished inner surface to minimize PE wear
- The alignment peg guides the metal liner during insertion, ensuring the metal liner is aligned with the acetabular cup

▲ PE Liner

- -Range of motion up to 200 degrees
- -The entry chamfer avoids impingement with the stem neck
- -The retention liner keeps the femoral head in place



SAFETY[™] **Bipolar Head**



- ▲ High polished Co-Cr-Mo Cup and Germany UHMWPE liner
- ▲ The special anti-dislocation snap ring is effective to prevent from dislocation
- ▲ The various sizes (with 1 mm increment) will meet different requirements

Cemented Acetabular Cup

- ▲ Horizontally & Vertically grooved backside provides rotational stability
- ▲ 22° Anti-dislocation Lip provides greater range of motion and minimizes dislocation

XLWzero Ceramic Liner



▲ Made of high-purity alumina-based and zirconiabased composite ceramic materials, providing a smoother surface finish, enhanced fracture resistance, and superior biocompatibility



Femoral Head

- ▲ Complete options in size: φ22, 24, 28, 32mm
- ▲ 2 options for material: stainless steel and Co-Cr-Mo







Ceramic Head







MYSMON™

Zirconium-niobium Femoral Head



- Excellent biocompatibilty
- Low MRI artifacts
- Fragmentation resistance
- Ceramic-like wear resistance
- Anti-fretting performance








Customization process

Consultancy

Image data collection \rightarrow 3D reconstruction of anatomical structure \rightarrow Medicine-engineering Interaction





Design

Product solutions → Prosthesis simulation assembly → Medicine-engineering Interaction







Simulation

3D printing bone model \rightarrow Medicine-engineering Interaction \rightarrow preoperative simulation operation with model





Customization

Image data collection \to 3D reconstruction of anatomical structure \to 3D printed bone model \to Surgical plan simulation \to Customized product design confirmation \to Submit for record \to 3D printed prosthesis or tools \to Delivery and surgery implementation \to Evaluation and follow-up







Hip Instrumentation


Total/Semi Hip Arthroplasty_Universal Instrumentation





3D SEE Instrument Set















DAA Special Tools (1)




Revision Stem Removal Instrumentation (I)



Revision Cup Removal Instrumentation



DAA Special Tools (2)



Revision Stem Removal Instrumentation (II)





JUST MEDICAL DEVICES (TIANJIN) CO., LTD

Add: No.27, Ziyang Road, Nankai District, Tianjin, China Tel: +86 22 2339 9501 Mobile: +86 178 2201 3519 (Whatsap

Web: www.just-ortho.com E-mail: sales@justmedical.cn

Version: 202407-15