

## **Shenzhen JCR Medical Technology Limited Company**



# **Customer First** Quality First





## Shenzhen JCR Medical Technology Limited Company

Add: Plant B, No.1, Tianfu Road, Tianliao Community, Yutang Street, Guangming District, Shenzhen, Guangdong, China. Fax:+86-755-82535606

WEB www.jcr-medibp.com jcrkj.en.alibaba.com jcrmedical.en.made-in-china.com

Mobile/Whatsapp/Wechat +86-13197290865 Email:sales01@jcr-med.cn/ Nancy.cao@jcr-med.com

Mobile/Whatsapp/Wechat +86-13471993401 Email:sales028@jcr-med.cn

#### Candy Tang

Mobile/Wechat +86-13926569635 Email:Candy.Tang@jcr-med.com

#### Christina Zou

Mobile/Whatsapp/Wechat +86-18071931385 Email:sales018@jcr-med.cn

#### Adele Xie

Mobile/Whatsapp/Wechat +86-15778609939 Email:sales06@jcr-med.cn



## Content

ECG Cables Series
Tourniquet Series
Water Trap
NIBP CUFF Series
Triple three-way valve device
Angiography Syringe
Surgical kits
Central Venous Catheter Kit
Pressure pump
Y-Valve



## **Company Profile**

Shenzhen JCR Medical Technology Limited Company founded in 2009, is located in the China Shenzhen special economic zone. We are manufacturing IBP transducers, SPO2 probes, Temperature probes, Tourniquets, Disposable Pressure Pump, ECG Cables, other medical tube and related accessories, related extension cables and related components. We are also distributor of the TE sensor.

Our products include IBP Transducers Series and Disposable Temperature Sensor Probes Series, Reusable Temperature Sensor Probes Series, Disposable SPO2 Sensor Probes Series, Reusable SPO2 Sensor Probes Series, compatible with all major brands of monitors IBP Cables Series, Temperature Extension Cables Series, SPO2 Extension Cables Series, Tourniquet, ECG Cables, NIBP Cuff, Pressure Pump, Medical tubing, Central Venous Catheter Kit Disposable Brain Electrode, Y-connector, manifolds, ECG electrode, Umbilical cord clamp, Drainage catheter kit, Ureteral stent, Heparin cap, Disposable cervical ripening ballon, catheter sheath, Suction tube and etc. We can provide customers with fully customized products, also provide our standard products.

Our factory covers an area of 12000 square meters, with 5000 square meters in accordance with China GMP standard 100000 grade clean workshop and 10000 grade clean laboratory. Our company is equipped with UV oven curing, accurate test system, super stable sealing machine and other equipment to produce finished products. At present, it can manufacture more than 400,000 sets per month. Mature production technology, good quality, and short lead time, 100% final test, to provide customers with ODM, OEM and customization services.









Our factory has high precision electronic horizontal injection molding machine and vertical injection molding machine, automatic packaging machine, blister packing machine and EO sterilization equipment, complete production line dedicated to sterile medical devices. Our company was rated as a national high-tech enterprise and technologically advanced 'little giant' enterprises in Guangdong Province. We also passed TUV ISO13485, CFDA and CE.

Our company adheres to the business philosophy of "people-oriented, technology development", adheres to the business policy of "customer first, quality first", focuses on the application of sensors in electronic medical devices, and is committed to providing customers with the best solutions. JCR Medical will always sincerely cooperate with domestic and foreign manufacturers as before, and look forward to working hand in hand with friends to make joint efforts for the development of the global medical industry.



### **ECG Cables Series**

#### Features

- ⊚3LD, 5LD, AHA, IEC, Clip, Snap
- © 5LD: AHA ( (white, green, brown, red, black), IEC (red, black, white, green, yellow)
- 3LD: AHA ( (white, red, black), IEC (red, green, yellow)
- © Each brand can do 5LD or 3 LD, 5 LD more than two buckle can measure breath.
- © Compatible monitors: Mindray MEC1000, MEC2000, PM7000, PM8000, PM9000, Bangjian, Edan, Goldway, Guoteng, Hexin Zhongdian, Genareal, Osen, Comen, Biolight, etc.

#### Disposable SpO2 Probe Product Reference Table: Food and Drug Administration 20210271

Connect patient the monitor cable connector ( Various types of monitor end connectors can be customized )		patient side connector				
		3LD Clip	5LD Clip	3LD Snap	5LD Snap	
Mindray 12P						
Nihon Kohden						
GE 11P				CHILL		
Comen 12P				1111	***************************************	
Philips 12P		- ann		100	1111	
Siemens 7P						
6P IEC						



### The application scenario of ECG lead wire products:



A—Connect the cable end of the monitor:

B—Main line cable, length. meter, material is TPU, total length 3.4 m:

**C**—HUB: **D**—patient cable, 0.9 meters: **E**—patient side connector:

## Complementary products for ECG patient cable ——Limb clips/chest suction balls

#### Limb clips

Multifunctional limb clip, Nickel plating, 4PCS/package, suitable for 3.0mm pin, 4.0mm banana plug Connect with the lead wire of 4.0mm snap button, compatible with all domestic and foreign ECG machines recommended compatible brands

#### chest suction ball

Multifunctional limb clip,Nickel plating,4PCS/package,suitable for3.0mm pin,4.0mm banana plug Connect with the lead wire of 4.0mm snap button, compatible with all domestic and foreign ECG machines recommended compatible brands





### **Tourniquet Series**

#### Scope of application

Mainly used for compression bandaging to eliminate lacunae and temporarily stop bleeding (non-arterial hemostasis).

#### **Features**

The ergonomic design ensures the normal blood flow of the ulnar artery and is convenient for fixation. Elastomer is soft, transparent and easy to observe Safe and comfortable; the elastic pressure belt is adjustable and easy to operate.







Snap Jiont type

Pneumatic tourniquet

### water trap

Rotate type



JJSP-01 Compatible with Mindray Dryline I Water Trap





### Blood pressure cuff series

#### reusable NIBP CUFF Series

#### **Product Specifications**

- Monitor NIBP CUFF series
- Reusable for different patients available
- $\bigcirc$  PU, nylon material with TPU , nylon are available , comfortable to touch and easy to clean
- © Fitting different brands of monitoring systems with variety of connectors equipped
- O If order large quantity we can OEM LOGO and print for free.









Brown PU leather NIBP cuff

Blue PU leather NIBP cuff

Black composite nylon material with TPU NIBP cuff Green gentle nylon NIBP cuff

## **Disposable NIBP Cuffs**

#### **Features**

- Monitor series NIBP Cuffs
- O Disposable/single patient use
- O Non-woven and TPU material for maximum patient comfort
- © Fitting different brands of monitoring series with variety of connectors equipped
- O If order large quantity we can OEM LOGO and print for free.



Disposable Non-woven NIBP Cuff



Disposable TPU NIBP Cuff, Neonate Series



## Triple three-way valve device

#### Introduction

Used for multi-channel infusion in various surgeries.



Smooth rotation, rotation, convenient multi-channel infusion.



High pressure three-way valve



Duplex three-way connection pipe



Triple three-way valve device

## Disposable Angiography Syringe

#### Introduction

In interventional procedures, it is used to precisely control the amount of injected contrast agent

#### **Features**



Disposable Angiography Syringe



Disposable Angiography Syringe Light purple

### Accurately control the injection dose, is easy to operate, and is easy to use.

## **Disposable Interventional Surgery Kit**

#### Introduction

With pre-installed extension lines and infusion device, help saving time. Luer lock design for secure connection Latex free



Disposable Interventional Surgery Kit



### **Disposable Central Venous Catheter Kit**

#### Scope of use

Intravenous infusion, manometry, and nutritional therapy were performed via internal jugular, subclavian, and femoral vein puncture.

#### **Product Features**

- ©Tube with good biocompatibility decreases infection rate, smooth inner and outer wall decreases thrombus formation and placement resistance, soft tip reduces injury.
- Connector with special structure is durable and anti-splitting; transparent extension tube.
- Tube is easy to observe, good elasticity makes it possible be opened and closed repeatedly, Y type needle and blue needle are optional.

#### **Model Specifications**

Туре	Catheter OD	Catheter Length(mm)	Guide wire (In*cm)	Populations commonly used
	14G	160、200、300	0.035*60	adult
Single	16G	160、200、300	0.035*60	adult
lumen	18G	130、160、200	0.021*60	pediatric
	20G	130、160、200	0.021*60	pediatric
Double	5Fr	130、160、200	0.021*60	pediatric
lumen	7Fr	160、200、300	0.035*60	adult
Triple	5.5Fr	130、160、200	0.021*60	pediatric
lumen	7Fr	160、200、300	0.035*60	adult



Single lumen Central Venous Catheter Kit

Single lumen Central Venous

Catheter Kit Developed View



Double lumen Central Venous Catheter Kit





Double lumen Central Venous Catheter Kit Developed View



Triple lumen Central Venous Catheter Kit Developed View



## Disposable Medical Balloon Pressure Pump

#### Application

For the purpose of dilating blood vessels or indwelling stents by pressurizing the balloon dilator catheter during PTCA surgery.

#### Feature

The maximum output pressure is 30ATM, the maximum capacity is 30ml for medical staff to use; Effective ergonomic design easy to operate, pressing type can be operated with one hand; The safety locking device can control the pressure and release the pressure quickly. The advance speed is even and the press is accurate.



JYLB-01 Disposable Medical Balloon Pressure Pump

## Type Y Connector

#### Introduction

Used for connecting guide tube in clinical operation and assisting in PTCA.

#### Feature

- © Rotating type: pressure-resistant and seal-fit operation. More convenient, tight safety, flexible valves.
- © Push Click type: Can be operated by one hand, with pressure tight sealing valve for rapid hemostasis, simple, flexible and safe operation, side wall with extension tube, more convenient operation.
- OSpring-type Y Hemostasis Valve



Push-type Y Hemostasis Valve

Push Click type Y Hemostasis Valve



Spring-type Y Hemostasis Valve



Y Hemostasis Valve Connector

T-type Y Hemostasis Valve



## Nasal Oxygen Cannula>>>



2.2m Transparent Nasal Oxygen Cannula

## Co2 Sampling Tube>>>



J00P0010 Sampling Tube

## Disposable body cavity tube>>>



body cavity tube

**Ureteral Stent** 

Pigtail Tube>>>



Ureteral Stent Pigtail Tube

## Spiral tube>>>



Spiral tube

## High Voltage Braid Tubing>>>



High Voltage Braid Tubing











